



Submission No.: PG01-9380

Session : Postgraduate Course 1 (Liver)

Date & Time, Place : November 16 (Thu), 08:30-10:00, Room 3F-1

Session Title : How do I do (RL donor): Hilar dissection & RL mobilization

Safe unilateral Hilar dissection for implementing new robotic donor programs

Yasir Alnema

King Faisal Specialist Hospital and Research Center, Saudi Arabia

Robotic approach provides higher magnification almost ten times. Fine articulation and wider range of motion provide accurate maneuvering and Safe dissection in narrowest space. This facilitates accurate dissection of the hepatic hilum, with indocyanine green (ICG) usage clear identification of the bile ducts therefore safe approach for bile duct division.

